


# CELT ACD

## A revolutionary approach to arterial closure

The Celt Arterial Closure Device (ACD) provides a unique and comfortable option for definitive closure using the existing procedure sheath.



Through the biocompatible 316 LVM implant, rapid haemostasis is achieved. Whilst also being clearly visible under fluoroscopy, safe and immediate re-access is possible if necessary.

- Definitive closure
- Rapid and simple deployment
- Uses existing sheath
- No skin incision
- No arterial dilation
- Visible under fluoroscopy making visualisation during deployment possible
- No patient discomfort
- No biological component

With the Celt ACD, you are 3 steps away from rapid haemostasis.


## 1, 2, 3 steps to arterial closure



1 Clockwise rotation of handle to open distal set of wings  
(inside puncture site)



2 Anti-clockwise rotation of handle to open proximal set of wings  
(outside of puncture site)



3 Release implant providing almost immediate haemostasis



### ORDERING DETAILS:

Product code	French Size	Box qty
KCLT-05	5	5
KCLT-06	6	5



**A faster, easier and more effective solution for your closure needs.**